

**FDA MEDICAL EQUIPMENT RECALLS AND ALERTS.** The following recalls are reported in Accordance with AFMLL 2-95, page CE-4, and should be treated as directed modifications in accordance with that AFMLL and AFI 41-201, Chapter 2, Section D. If you possess any of these items and have not received recall notification, you should contact the manufacturer for recall instructions. Suspension is only required for Class I items.  
**(FOM-P, Capt Paul J. Toth, DSN 343-7445)**

**CLASS I RECALLS:** None

**CLASS II RECALLS:**

**6515 NS**

**MDC 10134, 10144, 10145,  
97001**

**PRODUCT**

**CODE**

**MANUFACTURER**

**RECALLED BY**

**DISRIBUTION**

**QUANTITY**

**REASON**

**Anesthesia Units (General)**

Datex-Ohmeda AS/3 Anesthesia Delivery Unit (ADU), a complete anesthesia delivery unit with electronically controlled ventilator, gas delivery and agent vaporizing units. Recall #Z-004-1.

Serial Numbers: 40022203 and above.

Datex-OHMEDA, Inc., Tewksbury, Massachusetts.

Manufacturer, by voice mail on April 28, 2000, followed by visit beginning May 5, 2000. Firm-initiated field correction ongoing.

Nationwide and international.

1,150 units were distributed.

The device has reliability problems with fresh gas delivery, agent delivery and ventilation functions.

☐ None Present

☐ Action Taken \_\_\_\_\_

**6530 NS**

**MDC 16214**

**PRODUCT**

**CODE**

**MANUFACTURER**

**RECALLED BY**

**DISTRIBUTION**

**QUANTITY**

**REASON**

**Wheelchairs, Powered**

Power Chairs: a) Rascal Model Numbers RAS 545 and RAS 445;

b) Chauffeur Model Numbers: P48, S45, X48, and X48E.

Recall #Z-1029/1030-0.

Identified power chairs are: PC004492, PC004493, PC004494, PC004501, PC004502, PC004508, PC004511, PC004527, PC004532, PC004534, PC004536, PC004552, PC004553, PC004555, PC004556, PC004579, JS004709.

Electric Mobility Corporation, Sewell, New Jersey.

Manufacturer, by letter August 11, 2000. Firm-initiated recall ongoing.

a) North Carolina, Indiana, Florida, Pennsylvania, New York, Kentucky, California, Kansas, Mississippi, and West Virginia, England; b) Alabama.

17 chairs were distributed.

The Power chairs' controller program does not meet the original engineering specifications.

☐ None Present

☐ Action Taken \_\_\_\_\_

**6630 NS**  
**MDC 16405**  
PRODUCT

CODE  
MANUFACTURER  
SOFTWARE DEVELOPER  
RECALLED BY

DISTRIBUTION  
QUANTITY  
REASON

**Apheresis Unit**

Fenwal Autopheresis -C Plasmapheresis System, Models A-200, A-201, and A-401. Recall #B-860-0.

Product Code Numbers: 4R4550, 4R4560, 4R4561, R4R4585

Baxter Healthcare Corporation, Largo, Florida.

Baxter Healthcare Corporation, Round Lake, Illinois.

Baxter Healthcare Corporation, Deerfield, Illinois, by telephone and fax dated May 5, 2000, followed by letter dated June 15, 2000. Firm-initiated recall ongoing.

Nationwide and Canada.

5,250 units were distributed

The Autopheresis -C Plasmapheresis instruments with version 6.0 software may proceed to "Saline Rinse" without displaying the "Remove Plasma" prompt.

☐ None Present

☐ Action Taken \_\_\_\_\_

**CLASS II RECALLS:**

**6630 NS**  
**MDC 16817**  
PRODUCT

CODE  
MANUFACTURER

RECALLED BY

DISTRIBUTION

QUANTITY

REASON

**Blood Grouping Systems, Automated**

PK7100 and PK7200 CMV-PA Automated Blood Grouping Analyzer CMV Test System and Controls:

Olympus PK7100 Automated Blood Grouping System, Product Code 81KSZ;

Olympus PK7200 Automated Microplate System, Product Code 81KSZ;

Olympus PK CMV-PA Cytomegalovirus serological reagents, Product Code MZE.

Recall #B-1295-0.

All lot and serial numbers are affected.

PK7100 & PK7200 Automated Analyzer: Olympus Optical Co., Ltd., Tokyo,

Japan PK CMV-PA Systems and Controls: Fujirebio America, Inc., Fairfield, New Jersey; Fujirebio, Inc., Tokyo, Japan;

SCIMEDIX Corp., Denville, New Jersey.

Olympus America, Inc., Irving, Texas (distributor), by letter dated July 17, 2000. Firm-initiated recall ongoing.

North Carolina, California, Louisiana, Florida, Kentucky, Virginia, Illinois, Tennessee, Oklahoma, Canada.

PK7100: 2 analyzers; PK7200 38 analyzers were distributed. 412,000 tests per month were distributed.

Automated test systems, which may have produced incorrect CMV antibody test results.

☐ None Present

☐ Action Taken \_\_\_\_\_

**2. DRUG/SUPPLY PRODUCT RECALLS AND MEDICAL INFORMATION.** The Food and Drug Administration (FDA) advises that the drug products listed below were recalled by the manufacturers or distributors concerned. The FDA has classified these recalls as follows:

**CLASS I:** A situation in which there is a reasonable probability that the use of, or exposure to, a violative product will cause serious, adverse health consequences or death.

**CLASS II:** A situation in which the use of, or exposure to, a violative product may cause temporary or medically reversible adverse health consequences or where the probability of serious adverse health consequences is remote.

**CLASS III:** A situation in which the use of, or exposure to, a violative product is not likely to cause adverse health consequences.

Most of these items are nonstandard, and the possibility exists that medical activities may have purchased them locally. There is also a possibility that some of these items have NSNs assigned and may have been reported in earlier Q.A. messages. Activities having quantities of these items on hand will immediately suspend the materiel from issue and use. **CONUS** activities will contact the nearest office of the respective manufacturer or distributor for disposition instructions. **OVERSEAS** activities will report quantities suspended to AFMLO/FOM-P no later than **05JAN 01** for disposition instructions. They should include nomenclature; lot number; manufacturer; quantity suspended; strength size; requisition; and DSCP purchase order number, contract number, and stock record account number (SRAN). (FOM-P), **Bonnie Phillips DSN (343-4170)**

**CLASS I RECALLS: None**

**CLASS II RECALLS:**

NSN	6505 Nonstandard
PRODUCT	Source Leukocytes. Recall #B-1085-0.
CODE	Unit #3600010.
MANUFACTURER	Gulf Coast Regional Blood Center, Houston, Texas.
RECALLED BY	Manufacturer, by fax on May 5, 1999, and August 22 and 30, 1999. Firm-initiated recall ongoing.
DISTRIBUTION	Pennsylvania.
QUANTITY	1 unit was distributed.
REASON	Blood product was collected from a donor who had previously been deferred for high risk behavior.

[ ] None Present

[ ] Action Taken \_\_\_\_\_

NSN	6515 Nonstandard
PRODUCT	Centrex Bipolar Endoprosthesis Comp, used in conjunction with a femoral stem component in the reconstruction of a femoral head damaged by fresh fracture non-union aseptic necrosis of the femoral head and neck or by osteoarthritis or post-traumatic arthritis.
	Recall #Z-963-0.
CODE	Catalog No. 6225-2-846, Lot Code GOUDA.
MANUFACTURER	Howmedica Osteonics Corporation, Rutherford, New Jersey.
RECALLED BY	Manufacturer, by letter dated June 22, 2000. Firm-initiated recall ongoing.

DISTRIBUTION  
QUANTITY  
REASON

Nationwide, Australia and Canada.  
11 units were distributed.  
A 26mm UHWPE insert was issued instead of a 28mm insert.

☐ None Present  
☐ Action Taken \_\_\_\_\_  
\_\_\_\_\_

NSN  
PRODUCT

6550 Nonstandard  
Isoniazid drug lot No. 9335260 used in the following kits:  
BACTEC(tm) S.I.R.E. Drug Kit used for susceptibility testing of  
Mycobacterium tuberculosis, Catalog No. 442102, Lot Nos. 9327296,  
93442298, 9327298;  
BACTEC(tm) Isoniazid Drug Kit used for isoniazid susceptibility testing,  
Catalog No. 442146, Lot No. 9327297. Recall #Z-1032-0.

CODE  
MANUFACTURER  
RECALLED BY

See above.  
Becton Dickinson Microbiology Systems, Cockeysville, Maryland.  
Becton Dickinson Microbiology Systems, Sparks, Maryland, by letter on  
August 22 and 25, 2000, followed by telephone on September 6, 2000.  
Firm-initiated recall ongoing.

DISTRIBUTION

Nationwide, Hong Kong, Philippines, Australia, Brazil, New Delhi,  
Singapore, Puerto Rico, Belgium, Canada, Argentina.

QUANTITY  
REASON

148 Isoniazid Drug Kits and 490 S.I.R.E. Kits were distributed.  
Some bottles of the products listed above are labeled Isoniazid but  
actually contain Streptomycin.

☐ None Present  
☐ Action Taken \_\_\_\_\_  
\_\_\_\_\_

**CLASS III RECALLS:**

NSN  
PRODUCT

6505 Nonstandard  
BC(r) brand Headache Powder (650 mg Aspirin, 195 mg Salicylamide, 33.3  
Caffeine), packaged in 2's NDC# 1015800908, packaged in 6's NDC#  
1015800910, packaged in 24's NDC# 1015800912, packaged in 50's  
NDC# 1015800916. Recall#D-457-0.

CODE

Lot Numbers:  
M000580, M000581, M000582, M000583, M000857, M000858, M000859,  
M000860,  
M000861, M000862, M000863, M000864, M000865, M000866, M001703,  
M001704,  
M001705, M001706, M001707, M001708, M001709, M001710, M001711,  
M001712,  
M002205, M002206, M002207, M002208, M002209, M002210, M002211,  
M002212,  
M002213, M002214, M002771, M002772, M002773, M002774, M002780,  
M002903,  
M002904, M002905, M002906, M002907, M002908, M002909, M002910,  
M002911,

MANUFACTURER RECALLED BY  DISTRIBUTION QUANTITY   REASON	M002912, M003265, M003390, M003391, M003392, M003393, M003394, M003395, M003396, M003397, M003398, M003420, M003428, M003429, M003430, M003431, M003432, M003433, M003434, M003435, M003436, M003437, M004001, M004049, M004050, M004051, M004052, M004053. Block Drug Company, Memphis, Tennessee. Block Drug Company, Jersey, City, New Jersey, by letter on May 30, 2000, and June 23, 2000. Firm-initiated recall ongoing. Nationwide. BC Powder 2"s ---316,714/60dz BC Powder 6's ----180,232/60dz BC Powder 24's ----1,806/3dz BC Powder 50's ----18,188/3dz were distributed. Discoloring of glassine envelopes holding product.
---	--

[ ] None Present  
 [ ] Action Taken \_\_\_\_\_  
 \_\_\_\_\_

NSN UPDATE	6505 Nonstandard Jordan's Epinephrine Injection, USP, 1:1000, 1mg/mL, in 1 mL ampuls, which appeared in the May 24, 2000 Enforcement Report should have been listed as follows: Class II -- Recall #D-349-0, Lot Numbers: 980308, 980401, 980501, 980601, 990205; and Class III -- Recall #D-350-0, Lot Numbers: Lot Numbers: 980402, 990307, 990308.
---------------	---

[ ] None Present  
 [ ] Action Taken \_\_\_\_\_  
 \_\_\_\_\_

NSN PRODUCT   CODE MANUFACTURER RECALLED BY  DISTRIBUTION  QUANTITY REASON	6515 Nonstandard Replace 4.3mm Conical Abutment, 4.3mm diameter x 4.5 mm height, Catalog #43503. An abutment is used either for single or multiple unit, screw-retained prosthetic restorations. Recall #Z-1025-0. Lot #308039 EXP 2/05. Nobel Biocare USA, Inc., Yorba Linda, California. Manufacturer, by letter on June 13, 2000, followed by telephone. Firm-initiated recall ongoing. Florida, Illinois, Louisiana, Michigan, Minnesota, North Carolina, New York, Ohio, Pennsylvania, Texas, Hong Kong, Israel, Mexico. 40 units were distributed. The conical abutment was not the size as labeled.
---	---

[ ] None Present  
 [ ] Action Taken \_\_\_\_\_  
 \_\_\_\_\_

NSN 6515 Nonstandard  
PRODUCT IMC (Intramobile Connector) Conversion Screw, 3.3mm diameter, Catalog #8198CA. An IMC cylinder is used for either single or multiple unit, screw-retained prosthetic restorations. Recall #Z-1026-0.  
CODE Lot #310052.  
MANUFACTURER Nobel Biocare USA, Inc., Yorba Linda, California.  
RECALLED BY Manufacturer, by letter on June 14, 2000. Firm-initiated recall ongoing.  
DISTRIBUTION Arizona, Connecticut, Delaware, Florida, Indiana, Massachusetts, Maryland, Michigan, New Jersey, New York, Ohio, Texas.  
QUANTITY 29 units were distributed.  
REASON The screw in the package is incorrect.

☐ None Present  
☐ Action Taken \_\_\_\_\_  
\_\_\_\_\_

NSN 6515 Nonstandard  
PRODUCT Replace 4.3mm Conical Abutment, 4.3mm diameter x 1.5mm height, Catalog #43500. An abutment is used either for single or multiple unit, screw-retained prosthetic restorations. Recall Z-1028-0.  
CODE Lot #308028 EXP 2/05.  
MANUFACTURER Nobel Biocare USA, Inc., Yorba Linda, California.  
RECALLED BY Manufacturer, by letter on June 12, 2000. Firm-initiated recall ongoing.  
DISTRIBUTION California, Florida, Idaho, Indiana, New York, Colombia, Hong Kong, Mexico, South Korea.  
QUANTITY 41 units were distributed.  
REASON The conical abutment was not the size as labeled.

☐ None Present  
☐ Action Taken \_\_\_\_\_  
\_\_\_\_\_

NSN 6550 Nonstandard  
PRODUCT Vitros Immunodiagnostic Products brand of:  
Prolactin Reagent Pack, 100 test units per pack  
Catalog No.184 9793; and associated Vitros Prolactin Calibrators, Catalog No. 111 3596, for in vitro diagnostic use only. Recall #Z-1023/1024-0.  
CODE Catalog No. 1849793, Lot No.90, EXP 09/15/2000;  
Catalog No. 1112595, Lot No.90, EXP 09/15/2000.  
MANUFACTURER Ortho-Clinical Diagnostics, Inc., Forest Farm Estate Whitchurch, Cardiff, UK.  
RECALLED BY Ortho-Clinical Diagnostics, Inc., Rochester, New York, by letter dated May 24, 2000. Firm-initiated recall ongoing.  
DISTRIBUTION Alabama, Arizona, Kansas, Missouri, Virginia, Washington state, Argentina, Australia, Brazil, Canada, Chile, Colombia, India, Japan, Mexico, Puerto Rico, Panama, Singapore, Venezuela, England, France, Germany, Italy, and Spain.

QUANTITY  
REASON

764 packs were distributed.  
Calibration failures may occur when using this lot of product due to lower than expected light unit output for the level 1 calibrator.

☐ None Present  
☐ Action Taken \_\_\_\_\_  
\_\_\_\_\_

NSN  
PRODUCT

6550 Nonstandard  
ProCide-D Plus Sterilizing & Disinfecting Solution.  
Recall #Z-1031-0.

CODE  
MANUFACTURER  
RECALLED BY

Lot #200020002 EXP 7/02  
Metrex Research Corporation, Parker, Colorado.  
Sybron Dental Specialties, Inc., Orange, California, letter on August 7, 2000. Firm-initiated recall ongoing.  
Alabama, Iowa, Indiana, New York, Oregon, Pennsylvania, West Virginia, Texas.

DISTRIBUTION

QUANTITY  
REASON

280 gallons were distributed.  
Part of the subject lot was mislabeled as Procide-D Plus, which contains 3.4% glutaraldehyde, rather Procide-D, which contains 2.5% glutaraldehyde.

☐ None Present  
☐ Action Taken \_\_\_\_\_  
\_\_\_\_\_